
PRINTED THIN FILM YAZAKI **YTC America Inc. AMMONIA GAS SENSOR**

Detecting Ammonia at ppm-Level and High Relative Humidity

According to the World Heath Organization 10% of the world population is affected by chronic kidney disease (CKD) and millions die of it annually. Elevated levels of ammonia (NH₃) and trimethylamine (TMA) can in breath can be an indication of CKD, which may be treatable at an early stage. In contrast to blood panels detecting ammonia by a breathalyzer is less invasive and self tests can be performed at home. Typical ammonia levels are in a range of 0.5 ppm to 10 ppm and need to be detected at a relative humidity of 90%. Ammonia breath sensors need to be inert to exhaled oxygen, nitrogen and carbon-dioxide.

> **Screening For Chronic Kidney Disease (CKD) By Detecting Ammonia in Human Breath**







HIGHLIGHTS

Breathalyzer with single use ammonia detector strips

According to the WHO 10% of the world's population suffer from CKD and millions die annually


> Ammonia from the bloodstream is passed to the lungs

> > Spray Coater


Elevated Ammonia levels in breath can be a noninvasive method for early **CKD** detection

Dip-coat PANI on PET film


- Dopant A for NH₃
- Dopant B for TMA



Print Electrodes







Ammonia and TMA Sensor Setup with **On-Board Temperature and Humidity Sensor**

conceptualized


Ammonia (NH_3) and Trimethylamine (TMA) sensitive Polyaniline (PANI) films with interdigitated electrodes on PET substrates demonstrated

Prototype sensor electronics and software developed


1 to 10 ppm +/- 0.5 ppm detection of Ammonia at 90% relative humidity

Typical sensor response time of 30 seconds

Gas Manifold mass-flowcontrollers for controlled amounts of NH_3 / TMA in humid and CO₂-rich air



Typical sensing traces at 90% humidity: Initial slope is proportional to NH_3 / TMA gas concentration



Sensor Selectivity to NH₃ vs TMA

High slectivity between NH₃ and TMA with distinct dopants

No sensitivity to oxygen, nitrogen, and carbon-dioxide

CONTACT INFO: INFO@YTCA.COM

2021 © Copyright by YTC America, Inc. (YTCA). Information in this document is subject to change without notice. No reproduction is permitted without written approval from YTCA